


SDO1NZA1-XXXXX Valve Certifications and Testing:

USP VI HDPE Wetted Materials (USP <87>/<88>)
Assembled and 100% Leak Tested in ISO Class 7 Environment
USP <788> lot tested
USP <85> lot tested
ADI/TSE Free
BPA Free
Rated for 50 kGy ionizing radiation exposure
Contact Equibar for Shelf Life and USP <665> Extractable/Leachable Profile Details

DRAWING REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
All	R0	Created Drawing	05/24/2023	CD

Tolerances: (Unless otherwise noted)	DRAWN		5/24/2023					
	CHECKED		5/24/2023					
	RH		5/24/2023					
	QA							
	.x +/-0.030 .xx +/-0.015 .xxx +/-0.005	MFG			TITLE SDO1NZA1-XXXXX External			
Angular: +/- 1°	BG		5/24/2023					
All dimensions in Inches.	APPROVED							
				SIZE		DWG NO		REV
				C		SDO1NZA1-XXXXX External		R0
				SCALE	3 : 1		SHEET 1	OF 1
This drawing is Confidential and proprietary to Equibar, LLC. It may not be distributed or used for any purpose without written permission. Pressure regulator designs are protected by US Patent 6,886,591 and pending US and foreign patents.								